

Punctum Plugs



post@meditools.no
Tlf: 40 17 40 16



Long-Term Plugs

SUPEREAGLE®

The Simple Solution for Dry Eye




	PART #	DIAMETER	LENGTH
• Small, Medium & Large	31302	0.4 - 0.6mm (small)	1.4mm
• Super Comfort	31312	0.6 - 0.8mm (medium)	1.7mm
• Super Retention	31322	0.8 and larger (large)	1.7mm

Supplied: 2 individually packaged sterile plugs with inserter/dilator per twin pack.

SUPERFLEX®

The Comprehensive Solution for Dry Eye



	PART #	DIAMETER	LENGTH
• Custom Fit	31732	0.3mm	1.1mm
• Easy Insertion	31742	0.4mm	1.1mm
• Greater Patient Comfort	31752	0.5mm	1.2mm
	31762	0.6mm	1.3mm
	31772	0.7mm	1.4mm
	31782	0.8mm	1.5mm
	31792	0.9mm	1.7mm
	31702	1.0mm	1.8mm
	31712	1.1mm	2.0mm
	31722	1.2mm	2.0mm
	31802	1.3mm	2.0mm

Supplied: 2 individually packaged sterile plugs with inserter/dilator per twin pack.

EAGLEPLUG®

The original Tapered-Shaft™ punctum plug



	PART #	DIAMETER	LENGTH
• Designed to stay snugly in the punctal opening	30342	0.4mm	1.1mm
	30352	0.5mm	1.2mm
	30362	0.6mm	1.3mm
	30372	0.7mm	1.4mm
	30382	0.8mm	1.5mm

Supplied: 2 individually packaged sterile plugs with inserter/dilator per twin pack.

EAGLE FLEXPLUG®

The original flexible Tapered-Shaft™ punctum plug



	PART #	DIAMETER	LENGTH
• Designed for a flexible custom fit	30742	0.4mm	1.1mm
	30752	0.5mm	1.2mm
	30762	0.6mm	1.3mm
	30772	0.7mm	1.4mm
	30782	0.8mm	1.5mm
	30792	0.9mm	1.6mm

Supplied: 2 individually packaged sterile plugs with inserter/dilator per twin pack.

PLUG1™

The One-Size Solution for Dry Eye

	PART #	DIAMETER	LENGTH
• One size	30012	0.5mm-0.8mm	1.8mm

- Easy insertion
- Fits 0.5mm - 0.8mm puncta
- Easily reversible

Supplied: 2 individually packaged sterile plugs with inserter/dilator per twin pack.



Flow Controllers

EAGLEPLUG® TEARFLOW™

Up to a 150% larger lumen than our regular Flow Controller



	PART #	OUTSIDE DIA.	INSIDE DIA.	LENGTH
<ul style="list-style-type: none"> Designed for increased tear flow when compared to our original Tapered-Shaft Flow Controller 	3005	0.5mm	0.2mm	1.2mm
	3006	0.6mm	0.3mm	1.3mm
	3007	0.7mm	0.3mm	1.4mm
	3008	0.8mm	0.3mm	1.5mm
	3009	0.9mm	0.4mm	1.6mm
	3010	1.0mm	0.5mm	1.8mm

Supplied: 2 individually packaged sterile plugs with inserter/dilator per twin pack.

TAPERED SHAFT™ FLOW CONTROLLER

When Total Punctal Occlusion is too much



	PART #	OUTSIDE DIA.	INSIDE DIA.	LENGTH
<ul style="list-style-type: none"> Designed to help regulate lacrimal flow 	30652	0.5mm	0.2mm	1.4mm
	30662	0.6mm	0.2mm	1.4mm
	30672	0.7mm	0.2mm	1.7mm
	30682	0.8mm	0.2mm	1.9mm

Supplied: 2 individually packaged sterile plugs with inserter/dilator per twin pack.

Temporary Plugs

DURAPLUG®

Extended Temporary Inserts
Effective
60-180 days



Ideal for:	PART #	DIAMETER	LENGTH
<ul style="list-style-type: none"> Treatment of post LASIK induced dry eye Seasonal dry eye Retention of ocular medications 	0076	0.2mm	2.0mm
	0077	0.3mm	2.0mm
	0078	0.4mm	2.0mm
	0079	0.5mm	2.0mm

Supplied: 2 sterile inserts per package; 10 packages per box. Total of 20 inserts.
*This product is not CE Marked

COLLAGEN

Effective 3-5 days; dissolves completely in 7-10 days



	PART #	DIAMETER	LENGTH
<ul style="list-style-type: none"> A cost effective way to evaluate the benefits of permanent punctal occlusion. 	0025	Insert 0.2mm	2.0mm
	0026	Insert 0.3mm	2.0mm
	0027	Insert 0.4mm	2.0mm
	0052	"Shorts" 0.2mm	1.6mm
	0053	"Shorts" 0.3mm	1.6mm
	0054	"Shorts" 0.4mm	1.6mm

Supplied: 6 inserts per sterile foam pack, 10 packs per box. Total of 60 inserts.
*This product is not CE Marked

Gauge System

CORONEO™ PUNCTAL GAUGE

A Complete Gauging/Dilating System In One Instrument

- Graduated demarcated tips on both ends
- Gauge dilates from 0.5mm to 1.0mm diameters
- Can be used with all punctum plugs

PART#
2040

Supplied: 1 gauge per package.



*Developed in conjunction with Minas Coroneo, M.D., NSW, Australia

Tear Test

COLOR BAR™


Schirmer Tear Test*

- Soft, blue dye line
- Tear flow scale
- Pre-printed for right or left eye

PART #
0039

Supplied: 2 standardized sterile strips per pouch. 50 pouches per box. Total of 100 strips.

*This product is not CE Marked



Other Products

FORCEPS



DESCRIPTION:
Gold Punctal Plug Forceps
0.2mm Wide, Inside Groove
Pointed Tips

PART #
K5-6218



Jewelers, Textured Inner Surfaces At Tip

K5-6575



Jewelers No. 7, Sharply Curved, Single-Use

11-2007

SuperFlex® is a registered trademark of Katena Products
SuperEagle® is a registered trademark of Katena Products
EaglePlug® is a registered trademark of Katena Products
EaglePlug® TearFlow™ is a registered trademark of Katena Products
SuperFlex® is a registered trademark of Katena Products
Eagle FlexPlug® is a registered trademark of Katena Products
DuraPlug® is a registered trademark of Katena Products

Order Online: corzaeye.com

Order by Phone: 800-225-1195

For more information, visit corzaeye.com



Corza Ophthalmology (Katena Products, Inc.), 6 Campus Drive, Suite 310, Parsippany, NJ 07054, USA

Phone: 973.989.1600, Toll free (USA only): 800.225.1195, Fax: 973.989.8175

All trademarks are owned by the company and its affiliates.

©2022 Corza Medical

All Rights Reserved. REV-CAT-041117-Rev1